DE	H.E.L.P. Heparin 400,000 IU Gebrauchsanweisung – Bitte aufmerksam lesen!	2
GB	H.E.L.P. Heparin 400,000 IU Instruction for use – Please read carefully!	3
BG	H.E.L.P. Хепарин 400,000 IU Инструкции за ползване – моля, прочетете внимателно!	4
FR	Héparine H.E.L.P. 400,000 UI Mode d'emploi – À lire attentivement!	5
GR	Ηπαρίνη H.E.L.P 400,000 IU Οδηγίες Χρήσης – Παρακαλούμε διαβάστε προσεκτικά!	6
HR	H.E.L.P. Heparin 400,000 IU Upute za uporabu – pažljivo pročitati!	7
HU	H.E.L.P. Heparin 400,000 IU Alkalmazási útmutató – kérjük, olvassa el figyelmesen!	8
IT	Eparina H.E.L.P. 400,000 UI Istruzioni per l'uso – Leggere con attenzione!	9
LT	H.E.L.P. Heparinas 400,000 IU Naudojimo nuorodos – būtinai atidžiai perskaitykite!	10
PL	Heparyna H.E.L.P. 400,000 IU Instrukcja obsługi – proszę przeczytać uważnie!	11
RU	Гепарин H.E.L.P. 400,000 IU [ME] Инструкция по применению – прочтите внимательно!	12
SRB	H.E.L.P. Heparin 400,000 IU Uputstvo za upotrebu – molimo Vas da pažljivo pročitate!	13
TR	H.E.L.P. Heparin 400,000 IU Kullanım talimatı – Lütfen dikkatli okuyunuz!	14
CA/EN	H.E.L.P. Heparin 400,000 IU Instruction for use – Please read carefully!	15
CA/FR	Héparine H.E.L.P. 400,000 UI Mode d'emploi – À lire attentivement!	16

**B | BRAUN**

B. Braun Avitum AG
34209 Melsungen
Germany

123/470245/0615
REVISION DATE: 03/2015

470245/0615
432 x 296 mm
Plegado: 36 x 148 mm
Laetus: 1801
Black

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Lagerung bei

GB Instruction for use – Please read carefully!

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Pay attention to the indications and contraindications described in the instructions for use of the H.E.L.P. Futura set and in the operating manual of the Plasmat[®] Futura.

[4] Side effects

Pay attention to the side effects described in the instructions for use of the H.E.L.P. Futura set and in the operating manual of the Plasmat[®] Futura.

In rare cases, benzyl alcohol can cause hypersensitivity reactions.

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[6] Instructions for use

Unless prescribed otherwise, inject the contents of an H.E.L.P. Heparin 10,000 IU/mL vial (40 mL) into a bag of H.E.L.P. Acetate Buffer pH 4.85 (4000 ml) through the injection port immediately before use. Carefully mix the H.E.L.P. heparin sodium solution with the acetate buffer. Place the provided label "H.E.L.P. Heparin 400,000 IU " on the acetate buffer bag to show that heparin has been added. Connect the H.E.L.P. Acetate Buffer pH 4.85 solution to the H.E.L.P. Futura Set in accordance with the operating manual for the Plasmat[®] Futura.

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Comply with the specifications in the instructions for use of the H.E.L.P. Futura Treatment Set and in the operating manual of the Plasmat[®] Futura.

One vial (40 ml) of H.E.L.P. Heparin is typically used for each H.E.L.P. apheresis treatment.

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CA/EN Instruction for use – Please read carefully!

H.E.L.P. Heparin 400,000 IU

[1] Product description

Vial containing 40 ml of sterile heparin sodium solution for extracorporeal plasma treatment within the context of H.E.L.P. Apheresis. Only intended for extracorporeal use within the context of H.E.L.P. Apheresis. Sterile using aseptic processing techniques; endotoxin-free.

Composition:

1 ml contains:

Heparin sodium in accordance with the 4 th WHO standard (porcine intestinal mucosae)	10,000 IU
Sodium chloride	3 mg
Preservative: benzyl alcohol	10 mg
Water for injection	
Sodium hydroxide / hydrochloric acid for pH adjustment	

[2] Intended use

Sterile H.E.L.P. Heparin is part of the H.E.L.P. apheresis treatment set. It is intended to provide extracorporeal precipitation of LDL-cholesterol, lipoprotein(a), and fibrinogen within the context of H.E.L.P. therapy with the Plasmat[®] Futura.

[3] Indications and contraindications

Pay attention to the indications and contraindications described in the instructions for use of the H.E.L.P. Apheresis set and in the operating manual of the Plasmat[®] Futura Apheresis system.

[4] Side effects

Pay attention to the side effects described in the instructions for use of the H.E.L.P. Futura set and in the operating manual of the Plasmat[®] Futura Apheresis system. In rare cases, benzyl alcohol can cause hypersensitivity reactions.

[5] Warnings and precautions

- Not intended for intravenous infusion or subcutaneous use.
- H.E.L.P. Heparin must not be used after the expiration date stated on the receptacle and outer packaging.

- Only use if the solution is clear and colorless and if the receptacle is undamaged.
- Only remove the protective cap immediately before use.
- Any remaining solution must be discarded.
- Store the solution out of children's reach.

[6] Instructions for use

Unless prescribed otherwise, inject the contents of an H.E.L.P. Heparin 10,000 IU/mL vial (40 mL) into a bag of H.E.L.P. Acetate Buffer pH 4.85 (4000 mL) through the injection port immediately before use. Carefully mix the H.E.L.P. heparin sodium before use. Carefully mix the H.E.L.P. heparin sodium solution with the acetate buffer. Place the provided label "H.E.L.P. Heparin 400,000 IU" on the acetate buffer bag to show that heparin has been added. Connect the H.E.L.P. Acetate Buffer 4000 pH 4.85 solution to the H.E.L.P. Futura Apheresis system Set in accordance with the operating manual for the Plasmat[®] Futura Apheresis system.

[7] Method of administration and duration of treatment

Comply with the specifications in the instructions for use of the H.E.L.P. Apheresis Treatment Set and in the operating manual of the Plasmat[®] Futura Apheresis system. One vial (40 ml) of H.E.L.P. Heparin is typically used for each H.E.L.P. apheresis treatment.

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