

Askina® DresSil

Product Data Sheet



Administrative information

Legal Manufacturer Winner Medical Co. Ltd. Winner Industrial Park, No 660 Bulong Road, Longhua District, 518109 Shenzhen, China

Importer B. Braun Avitum AG, Schwarzenberger Weg 73-79, 34212 Melsungen, Germany

Product management Wound Management

Last Update 30/05/2025

Description, composition and properties of the device

Trade Name Askina® DresSil

Reference WIN5291010; WIN5291505; WIN5292005

Medical Class Class IIb (MDR)

Description of the device Askina DresSil®, a silicone foam dressing without border, combines the superior absorption of foam and the gentle adhesion of silicone. It absorbs exudate through silicone perforations from the wound and provides a moist wound healing environment to promote wound healing. The silicone wound contact layer reduces pain and trauma during dressing changes and enhances patient's comfort. The waterproof outer layer protects the wound from dirt and bacteria.

Composition of the device Polyurethane film backing
Polyurethane foam
Silicone net
Release film

Key drivers & Indications	
Key drivers	Maintains a moist wound environment Minimizes pain during dressing changes Easy application and removal Breathable Does not adhere to the wound bed
Intended Purpose	Askina DresSil® is a highly conformable dressing that absorbs exudate and maintains a moist wound environment to promote wound healing.
Indications	Askina DresSil® is designed for a wide range of exuding wounds such as leg and foot ulcers, pressure ulcers and traumatic wounds, e.g. skin tears and secondary healing wounds.
Precautions of use	<p>Askina® DresSil is contraindicated in:</p> <ul style="list-style-type: none">• Check the wound for signs of infection before use, if infection occurs, see a health care professional.• Do not use on the patients with a known hypersensitivity to the product itself or to its components.• Do not use on third degree burns. <p>Warnings:</p> <ul style="list-style-type: none">• Do not reuse. Reuse will cause cross-contamination.• Do not use if package is damaged or open.• Do not use together with oxidizing agents such as hypochlorite solutions or hydrogen peroxide.• Due to the gentle adhesive used, do not use this dressing to secure other medical devices to the patient.• If reddening or sensitization occurs discontinue use and see a healthcare professional.
IFU: Yes/No	Yes
Reusable/single use device	Single use device
Sterilization process	
Sterile: Yes/No	Yes
Sterilization method	Ethylene Oxide (EO)
Conservation and storage conditions	
Storage conditions	The product should be stored in dry conditions. Keep the product away from direct sunlight and keep dry.
Transport conditions	The product should be transported in dry conditions. Keep the product away from direct sunlight and keep dry.
Shelf life	3 years

Safety in use

Technical: MRI, X-ray detectable	N/A
Biocompatibility	Biocompatibility studies showed the device has no indication of eliciting the following responses: Cytotoxicity Irritation or sensitization Systemic Toxicity Local effects after implantation

Standards List

EN ISO 13485:2016	Medical devices – Quality management systems – Requirements for regulatory purposes
EN ISO 14971: 2019	Medical devices – Application of risk management to medical devices
EN ISO 20417:2021	Information supplied by the manufacturer of medical devices
EN ISO 15223-1: 2021	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
ISO 15223-2: 2010	Medical devices – Symbols to be used with medical device labels, labelling, and information to be supplied – Part 2: Symbol development, selection and validation
EN ISO 10993-1: 2020	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process
EN ISO 10993-5: 2009	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-10: 2013	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization
EN ISO 10993-11: 2018	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
EN ISO 10993-6: 2016	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
EN ISO 10993-11: 2018	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
EN ISO 11135: 2014/ A1: 2019	Sterilization of health-care products – Ethylene oxide – Requirements for the development, validation and routine control of a sterilization process for medical devices – Amendment 1: Revision of Annex E, Single batch release
ISO 10993-7: 2008/Amd 1: 2019	Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals – Amendment 1: Applicability of allowable limits for neonates and infants

Packaging & references			
Reference	Description	Size	Box Quantity
WIN5291010	Askina® DresSil	10 x 10 cm (non-border)	10
WIN5291505	Askina® DresSil	15 x 15 cm (non-border)	5
WIN5292005	Askina® DresSil	20 x 20 cm (non-border)	5