

Actreen® Lite Cath

Simplify every day life

NEW



Continence Care

Self intermittent catheterisation

B | BRAUN
SHARING EXPERTISE

Simplify every day life

Ready to use concept

User friendly design

- No additional wrapping
- Handy packaging
 - Two ways of opening the box
 - Easy accessibility to the product

- Easy handling
 - Two finger grab eyelets enabling opening with one finger only


- Universal connector
 - with an enlarged base for an easy grip

No Touch and Ready to Use


- No touch system
 - Maintained sterility during catheterisation
 - No direct contact between fingers and catheter
- Discreet and sterile catheterisation, wherever you are
 - High resistant film
 - Additional precut opening
 - Simplifies access to catheter

Men: one 1/3 part from the catheter tip
→ Release first part of the catheter from protective sheath

Women: close to the connector
→ Release the catheter from protective sheath completely



Actreen®, a catheter developed by patient for patient






Actreen® Lite Cath



Pack of
30 sterile catheters



		CH	REFERENCE
Actreen® Lite Tiemann 45 cm		CH08	228108
		CH10	228110
		CH12	228112
		CH14	228114
		CH16	228116
		CH18	228118
Actreen® Lite Nelaton 45 cm		CH08	228208
		CH10	228210
		CH12	228212
		CH14	228214
		CH16	228216
Actreen® Lite Nelaton 20 cm		CH06	228306
		CH08	228308
		CH10	228310
		CH12	228312
		CH14	228314
		CH16	228316