

Urimed[®] Bag

The proof of reliable quality



Adult urine collection bags

Urimed® Bag Range

The Urimed® Bag range meets patient's needs to have a **reliable** product adapted to each application.

We provide patients with, **high quality** urine collection bags which are **easy to use** and **comfortable** for their daily activities.

COMPLETE Product Range

B. Braun offers a complete range of adult urine collection bags including leg and bedside bags.

The range includes:

- Urimed® leg bags
- Urimed® night & day bags
- Urimed® sampling port
- A complementary range of Urimed® accessories

HIGH QUALITY Bags

Our wide range of urine bags meet European standard requirements.

- Urimed® leg bags have a soft back foil, avoiding skin irritation for better comfort.
- The Urimed® Bag products have a non kinking tube and large diameter for a better flow of urine.



Free Flow Lever Tap



Cross Outlet with stopper


SECURE Drainage

The drainage outlets are designed to be easily used for all levels of hand dexterity.

- Urimed® leg bags have a free flow lever tap for an easy and secure drainage with one hand.
- Urimed® night & day bags and sampling port have a cross outlet with stopper.

ANTI REFLUX
VALVE
ON ALL BAGS

Urimed® leg bags



● **Urimed® Bag Plus 500 ml**

- Non sterile
- 60 cm length tube, cut to fit

For a connection to
**B. Braun MEC*
Urimed® Vision**



*Male external catheter


● **Urimed® Tribag Plus 800 ml**

- Specific 3 chamber pouch
- Non sterile and sterile
- 60 cm length tube




● **Urimed® Bag Plus 1,5L**

- Non sterile
- 90 cm tube, cut to fit



For a connection to
**B. Braun urostomy
pouch:
Almarys® twin + Uro**



Urimed® night & day bags



• Urimed® Bag 2L

- Non sterile
- 120 cm tube




COST
EFFECTIVE

• Medicare 2L

- Sterile and non sterile
- 90 cm length tube

Urimed® sampling port



AIR FILTER FOR
A COMPLETE FILLING
OF THE BAG

New

• Urimed® SP 2L

- Sterile
- 100 cm length tube
- Clamp
- Urine sampling with **needle free** syringe
→ Eliminate the risk of needle stick injury



Hygienic foldable drain outlet




Needle free sample port

Urimed® Accessories

LEG FIXATION


Urimed® Fix

- Can be worn on the thigh or lower leg




Urimed® Straps

- Adapted for leg bag



Urimed® Klett

- Fixation for tube and catheter



BED FIXATION

Urimed® Handle

- To hang at the bedside or wheelchair



Leg Bags

	Length tube	Tube Ø (mm)	Box	References
Urimed® Bag Plus Non sterile	60 cm cut to fit	9 mm	10	28501*
Urimed® Tribag Plus 800 ml – 3-Chamber Sterile	60 cm connected connector	9 mm	10	28305*
Non sterile	60 cm cut to fit	9 mm	10	28306*
Urimed® Bag Plus 1,5L – Non sterile	90 cm cut to fit	9 mm	10	28150*

Night & Day Bags

	Length tube	Tube Ø (mm)	Box	References
Urimed® Bag 2L Non sterile	120 cm connected connector	9 mm	30	28300*
Medicare® 2L Sterile	90 cm	6,5 mm	10	28115*
Non sterile	90 cm	6,5 mm	10	28114*

Bags with sampling Port

	Length tube	Tube Ø (mm)	Box	References
Urimed® SP 2L	100 cm connected connector	9 mm	10 50	28610* 28650

Accessories

	Size	Units/box	References
Urimed® Fix**	small	2	68520*
	medium	2	68530*
	large	2	68540*
Urimed® Straps		1 pair	68550*
Urimed® Klett	small	10	28260*
	large	10	28261*
Urimed® Handle		1	68548*

* For each country, a specific letter is added at the end of the reference radical. ** Distributed by B. Braun