Urimed[®] SP Range

Quality within everyone's reach



Urine bags with needle free sample port





Urimed[®] SP

Urine bag with needle free sample port




Urimed[®] SP+

Urine bag with needle free sample port & hanging system integrated



Urimed[®] SP+



Urimed [®] SP+	Length tube	Tube diameter	REFERENCE	UNITS/BOX
O	130 cm	9 mm	28630	30

Developed and manufactured by: B.Braun Medical S.A.S. | RCS Nanterre 562050856 Urinary Care 204, avenue du Maréchal Juin | 92100 Boulogne-Billancourt Cedex France www.bbraun.com

The content of this document, in particular the data, the information, the trademarks and the logos are B BRAUN MEDICAL SAS's sole property and are protected under copyright 2014, B BRAUN MEDICAL, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of B BRAUN MEDICAL SAS's copyright and other proprietary rights.

The products commercialized by B BRAUN MEDICAL shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes.