



SPINE SURGERY

AESCULAP[®] CCR SYSTEM
TITANIUM
CASPAR[®] CERVICAL RETRACTOR SYSTEM

AESCULAP® CCR SYSTEM TITANIUM

FEATURES



EXTENSIVE BLADE RANGE

- 19 mm system for microsurgical techniques (green color-coded)
- 24 mm system for standard techniques (blue color-coded)
- blade lengths from 25 mm up to 75 mm in different tooth designs

SIDE-LOAD MECHANISM

- smart intra-operative handling
- modular handles to be used as hand-held retractors
- compatible with previous ball snap closure

DOUBLE HINGE RETRACTORS

- allows for intra-operative adaption
- basic retractor with tiltable butterfly screw

FENESTRATED TITANIUM BLADES

- enhanced grip and fit
- designed to allow for a better view and less weight
- semi-radiolucent
- color-coded according to blade width



AESCULAP® CCR SYSTEM TITANIUM

PRODUCT INFORMATION

CCR TITANIUM BLADES



toothed blades
(medial / lateral)

Art. No.	Size	Art. No.
BW050T	19 x 25 mm	BW060T
BW051T	19 x 30 mm	BW061T
BW052T	19 x 35 mm	BW062T
BW053T	19 x 40 mm	BW063T
BW054T	19 x 45 mm	BW064T
BW055T	19 x 50 mm	BW065T
BW056T	19 x 55 mm	BW066T
BW057T	19 x 60 mm	BW067T



blunt blades
(cranial-caudal)



lateral blades
(small teeth)

Art. No.	Size	Art. No.
BW141T	24 x 30 mm	BW151T
BW142T	24 x 35 mm	BW152T
BW143T	24 x 40 mm	BW153T
BW144T	24 x 45 mm	BW154T
BW145T	24 x 50 mm	BW155T
BW146T	24 x 55 mm	BW156T
BW147T	24 x 60 mm	BW157T
BW148T	24 x 65 mm	BW158T
BW149T	24 x 70 mm	BW159T



medial blades
(large teeth)



cranial-caudal blades
(blunt)

Art. No.	Size
BW162T	24 x 35 mm
BW163T	24 x 40 mm
BW164T	24 x 45 mm
BW165T	24 x 50 mm
BW166T	24 x 55 mm
BW167T	24 x 60 mm
BW168T	24 x 65 mm
BW169T	24 x 70 mm
BW170T	24 x 75 mm

RETRACTORS


BW048R

CCR counter retractor, preferably used for multi-level retraction



BW049R

CCR basic cervical retractor, preferably used for transverse retraction, double hinge, tiltable butterfly screw



BW047R

CCR basic counter retractor, preferably used for cranial-caudal retraction, double hinge



BW010R

Handle for blades to be used as hand-held retractor



BW011P

Depth gauge to determine the appropriate blade length



AESCU LAP[®] CCR SYSTEM TITANIUM

RETRACTOR SETS

BW124 – 24 mm TITANIUM BLADE SYSTEM

BW124 – containing		
Art. No.	Units	Description
BW125P	1	CCR tray, with lid
BW049R	1	CCR basic retractor
BW047R	1	CCR basic counter retractor
BW048R	1	CCR counter retractor
BW010R	2	Handle
BW011P	1	Depth gauge
BW141T	1	CCR blade, lateral 24 x 30 mm
BW142T	1	CCR blade, lateral 24 x 35 mm
BW143T	1	CCR blade, lateral 24 x 40 mm
BW144T	1	CCR blade, lateral 24 x 45 mm
BW145T	1	CCR blade, lateral 24 x 50 mm
BW146T	1	CCR blade, lateral 24 x 55 mm
BW147T	1	CCR blade, lateral 24 x 60 mm
BW148T	1	CCR blade, lateral 24 x 65 mm
BW149T	1	CCR blade, lateral 24 x 70 mm
BW151T	1	CCR blade, medial 24 x 30 mm
BW152T	1	CCR blade, medial 24 x 35 mm
BW153T	1	CCR blade, medial 24 x 40 mm
BW154T	1	CCR blade, medial 24 x 45 mm
BW155T	1	CCR blade, medial 24 x 50 mm
BW156T	1	CCR blade, medial 24 x 55 mm
BW157T	1	CCR blade, medial 24 x 60 mm
BW158T	1	CCR blade, medial 24 x 65 mm
BW159T	1	CCR blade, medial 24 x 70 mm
BW162T	1	CCR blade, blunt 24 x 35 mm
BW163T	1	CCR blade, blunt 24 x 40 mm
BW164T	1	CCR blade, blunt 24 x 45 mm
BW165T	1	CCR blade, blunt 24 x 50 mm
BW166T	1	CCR blade, blunt 24 x 55 mm
BW167T	1	CCR blade, blunt 24 x 60 mm
BW168T	1	CCR blade, blunt 24 x 65 mm
BW169T	1	CCR blade, blunt 24 x 70 mm
BW170T	1	CCR blade, blunt 24 x 75 mm



BW122 – Basic retractor set containing		
Art. No.	Units	Description
BW049R	1	CCR basic retractor
BW010R	2	Handle
BW143T	1	CCR blade, lateral 24 x 40 mm
BW144T	1	CCR blade, lateral 24 x 45 mm
BW145T	1	CCR blade, lateral 24 x 50 mm
BW146T	1	CCR blade, lateral 24 x 55 mm
BW147T	1	CCR blade, lateral 24 x 60 mm
BW153T	1	CCR blade, medial 24 x 40 mm
BW154T	1	CCR blade, medial 24 x 45 mm
BW155T	1	CCR blade, medial 24 x 50 mm
BW156T	1	CCR blade, medial 24 x 55 mm
BW157T	1	CCR blade, medial 24 x 60 mm

BW119 – 19 mm TITANIUM BLADE SYSTEM FOR MICROSURGICAL TECHNIQUES

BW119 – containing

Art. No.	Units	Description
BW125P	1	CCR tray, with lid
BW049R	1	CCR basic retractor
BW047R	1	CCR basic counter retractor
BW048R	1	CCR counter retractor
BW010R	2	Handle
BW011P	1	Depth gauge
BW050T	2	CCR blade, toothed 19 x 25 mm
BW051T	2	CCR blade, toothed 19 x 30 mm
BW052T	2	CCR blade, toothed 19 x 35 mm
BW053T	2	CCR blade, toothed 19 x 40 mm
BW054T	2	CCR blade, toothed 19 x 45 mm
BW055T	2	CCR blade, toothed 19 x 50 mm
BW056T	2	CCR blade, toothed 19 x 55 mm
BW057T	2	CCR blade, toothed 19 x 60 mm
BW060T	1	CCR blade, blunt 19 x 25 mm
BW061T	1	CCR blade, blunt 19 x 30 mm
BW062T	1	CCR blade, blunt 19 x 35 mm
BW063T	1	CCR blade, blunt 19 x 40 mm
BW064T	1	CCR blade, blunt 19 x 45 mm
BW065T	1	CCR blade, blunt 19 x 50 mm
BW066T	1	CCR blade, blunt 19 x 55 mm
BW067T	1	CCR blade, blunt 19 x 60 mm



BW123 – Basic counter retractor set containing

Art. No.	Units	Description
BW047R	1	CCR basic counter retractor
BW164T	1	CCR blade, lateral 24 x 45 mm
BW165T	1	CCR blade, lateral 24 x 50 mm
BW166T	1	CCR blade, lateral 24 x 55 mm
BW167T	1	CCR blade, lateral 24 x 60 mm
BW168T	1	CCR blade, lateral 24 x 65 mm

AESFULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "AESFULAP" and the product trademark "CASPAR" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.